CORRESPONDENCE

Inside Alzheimer brain with CLARITY: senile plaques, neurofibrillary tangles and axons in 3-D

Kunie Ando · Quentin Laborde · Adina Lazar · David Godefroy · Ihsen Youssef · Majid Amar · Amy Pooler · Marie-Claude Potier · Benoit Delatour · Charles Duyckaerts

Received: 21 June 2014 / Revised: 9 July 2014 / Accepted: 9 July 2014 / Published online: 29 July 2014 © The Author(s) 2014. This article is published with open access at Springerlink.com

Alzheimer's disease (AD) is characterized by extracellular deposits of amyloid β (A β) peptide and intracellular tau aggregates. Examination of these lesions in *post-mortem* brain tissue routinely provides a two-dimensional picture of the pathology with a limited depth of field. Several tissue-clearing protocols have been recently published [1, 4–6, 8, 10, 11, 16]. We compared the quality of immunostaining on *post-mortem* human brain tissues using several clearing techniques including CLARITY [4, 5], Scale [8], SeeDB [9], and 3DISCO [6]. Results were best with CLARITY. In that method, acrylamide hydrogel links the protein of the tissue and maintains its structure. Removing the lipids makes the block transparent and facilitates

K. Ando and Q. Laborde contributed equally.

Electronic supplementary material The online version of this article (doi:10.1007/s00401-014-1322-y) contains supplementary material, which is available to authorized users.

K. Ando · Q. Laborde · A. Lazar · C. Duyckaerts Laboratoire de Neuropathologie Escourolle, AP-HP, Hôpital de la Pitié Salpêtrière, 75013 Paris, France

K. Ando · Q. Laborde · A. Lazar · I. Youssef · M. Amar · M.-C. Potier · B. Delatour · C. Duyckaerts (⊠) Alzheimer's and Prion disease team, UPMC Université Paris 06 UMR S 1127, Inserm, U 1127, and CNRS UMR 7225, ICM, 75013 Paris, France e-mail: charles.duyckaerts@psl.aphp.fr

D. Godefroy

UMR S 968 Inserm/UPMC/CHNO des Quinze-Vingts, Centre de Recherche Institut de la Vision, 75012 Paris, France

A. Pooler

Department of Neuroscience, Institute of Psychiatry, King's College London, DeCrespigny Park, London SE5 8AF, UK

antibody diffusion [4, 5]. We noticed that the results were improved if ScaleA2 [8] was used to mount the block after CLARITY.

Formalin-fixed frontal cortex samples from two controls and five AD individuals (Braak VI and Thal 5) were obtained from the Brain Bank GIE NeuroCEB with legal consent. Brain samples were cut at a thickness of 500 µm with a vibratome and were processed by the CLARITY technique [5]. Embedded in an acrylamide hydrogel constituted of 4 % paraformaldehyde (PFA), 4 % acrylamide, 0.25 % temperature-triggering initiator VA-044, PBS, the tissues were passively clarified at 37 °C for 2 weeks in the clearing solution (200 mM boric acid, 4 % w/v SDS, pH 8.5). The blocks were immunostained with a rabbit polyclonal anti-tau B19 antibody [3] and a mouse monoclonal anti-Aß 4G8 antibody (Covance). Alexa Fluor® 488 goat anti-rabbit and Alexa Fluor® 568 goat anti-mouse antibodies (Life technologies) were used as secondary antibodies. For AB and neurofilament double immunohistochemistry, the tissues were first incubated with mouse monoclonal anti-neurofilament antibody (M0762 Dako) and then with Alexa Fluor[®] 568 goat anti-mouse antibody. After a blocking stage (10 % normal mouse serum), the tissues were incubated with biotin-labelled 4G8 anti-AB antibody (Covance), revealed with DyLight 488-labelled streptavidin (KPL, Eurobio, France). Triple staining (Aβ, tau and neurofilament) was also performed: the block was incubated with anti-tau B19 and anti-neurofilament M0762 antibodies revealed with the secondary Alexa Fluor® 488 goat antirabbit and Alexa Fluor[®] 568 goat anti-mouse antibodies. The third antibody was a biotin-labelled anti-AB 4G8 antibody (Covance) revealed with streptavidin-Alexa 405 (Life technology). The mounting protocol was modified from the original method: we added a fixation step (4 % PFA in PBS for 15 min) at the end of the procedure to improve the stability of the immunostaining, and an incubation step (0.2 M glycine for 15 min) to quench autofluorescence. We used ScaleA2 solution as mounting medium [8], instead of 80 % glycerol or FocusClear described in the original method [5]. ScaleA2 increased the transparency of the tissue compared to 80 % glycerol, lowered the cost of the experiment compared to FocusClear and shortened incubation time (overnight compared to 2 days). The immuno-labellings were analysed with an upright confocal microscope (Olympus Fluoview Fv1000). The Z-stack images were re-constructed using Imaris software (Bitplane).

The volume of the "clarified" block was increased by approximately 40 % after passive removal of the lipids. While 80 % glycerol or FocusClear restored the volume in the original CLARITY technique, ScaleA2 further increased it but kept the morphological quality of the technique. The combination of CLARITY and ScaleA2 is not appropriate for stereological measures because of the volume expansion. Despite the increase in volume, the axons (shown by neurofilament immunostaining) remained continuous (Supplementary videos 2-3). The clarified tissues were fully transparent. The video, obtained with confocal microscopy, provides a 3-D view of the lesions revealed by $A\beta$ and tau immunohistochemistry (video 1). Amyloid deposits appeared to have diverse 3-D structure in AD brains: some deposits were dense and focal; some diffuse deposits appeared hollow (video 1). Accumulation of AB was regularly spaced in clumps throughout the cortex (video 1). Tau accumulated in segmented and discontinuous neuritic processes. Many dystrophic neurites were associated with dense, focal, presumably "mature" A β focal deposits [2].

Plaque-induced distorted neurites have been previously reported by in vivo observation in AD transgenic animal [7, 12, 13, 15] and in the affected areas of *post-mortem* human AD brains [10, 14]. Our video of A β and neurofilament double immunostaining showed the persistence of a large number of axons, one of them (traced in yellow) deflected by the focal deposits (video 2).

Finally, a triple immunostaining for $A\beta$ (blue), tau (green) and neurofilament (red) showed that the tangle bearing neuron that we identified had no visible axonal projection (video 3). The region of the $A\beta$ focal deposit contained few axons. The $A\beta$ immunoreactivity (revealed with an avidin biotin amplification step) was excellent through the whole thickness of the block. Neurofilament and tau immunoreactivities (revealed with no amplification step) were stronger in a 150-µm-thick layer at the top of the block. The transparency of the block was, however, homogeneous. The 150-µm-thick layer of high immunoreactivity was detected at the top and at the bottom of the block: it was not related to a defect in the penetration of the laser beam during the imaging process but to better antibody penetration near the surface. Longer incubation

time, higher concentration of detergent in the buffer and various amplification systems will be tested in the future to improve the homogeneity of the signal throughout the thickness of the sample.

These first results uncover the 3-D architecture of AD lesions and demonstrate the utility of the CLARITY technique in studying the distribution of AD pathology in 3-D, even in *post-mortem* human tissue.

Acknowledgments The research leading to these results has received funding from the program "Investissements d'avenir" ANR-10-IAIHU-06. Supported by grants ANR Franco-Canada CholAD (U975-P-A1L4-RP-RPV10014DDA-R10178DD) and MassImage (U975-P-A1L4-RP-RPV10015DDA-R10179DD). The financial help of Association France Alzheimer is greatly acknowledged. KA is a grant recipient of French "Fondation Plan Alzheimer". Rabbit polyclonal anti-tau B19 antibody was a generous gift from Prof. Jean-Pierre Brion, Free University of Bruxelles, Belgium. The brain tissues were kindly provided by the GIE NeuroCEB brain bank (France Alzheimer, France Parkinson, ARSEP, CSC). Analysis by confocal microscopy was performed at the cellular imaging platform of Pitié-Salpêtrière PICPS.

Open Access This article is distributed under the terms of the Creative Commons Attribution License which permits any use, distribution, and reproduction in any medium, provided the original author(s) and the source are credited.

References

- Becker K, Jahrling N, Saghafi S, Weiler R, Dodt HU (2012) Chemical clearing and dehydration of GFP expressing mouse brains. PLoS One 7:e33916
- Braak E, Braak H, Mandelkow EM (1994) A sequence of cytoskeleton changes related to the formation of neurofibrillary tangles and neuropil threads. Acta Neuropathol 87:554–567
- Brion JP, Hanger DP, Bruce MT, Couck AM, Flament-Durand J, Anderton BH (1991) Tau in Alzheimer neurofibrillary tangles. Nand C-terminal regions are differentially associated with paired helical filaments and the location of a putative abnormal phosphorylation site. Biochem J 273(Pt 1):127–133
- Chung K, Deisseroth K (2013) CLARITY for mapping the nervous system. Nat Methods 10:508–513
- Chung K, Wallace J, Kim SY et al (2013) Structural and molecular interrogation of intact biological systems. Nature 497:332–337
- Erturk A, Becker K, Jahrling N et al (2012) Three-dimensional imaging of solvent-cleared organs using 3DISCO. Nat Protoc 7:1983–1995
- Garcia-Alloza M, Borrelli LA, Rozkalne A, Hyman BT, Bacskai BJ (2007) Curcumin labels amyloid pathology in vivo, disrupts existing plaques, and partially restores distorted neurites in an Alzheimer mouse model. J Neurochem 102:1095–1104
- Hama H, Kurokawa H, Kawano H et al (2011) Scale: a chemical approach for fluorescence imaging and reconstruction of transparent mouse brain. Nat Neurosci 14:1481–1488
- Ke MT, Fujimoto S, Imai T (2013) SeeDB: a simple and morphology-preserving optical clearing agent for neuronal circuit reconstruction. Nat Neurosci 16:1154–1161
- 10. Knowles RB, Wyart C, Buldyrev SV et al (1999) Plaqueinduced neurite abnormalities: implications for disruption of

neural networks in Alzheimer's disease. Proc Natl Acad Sci USA 96:5274–5279

- Kuwajima T, Sitko AA, Bhansali P, Jurgens C, Guido W, Mason C (2013) ClearT: a detergent- and solvent-free clearing method for neuronal and non-neuronal tissue. Development 140:1364–1368
- Le R, Cruz L, Urbanc B et al (2001) Plaque-induced abnormalities in neurite geometry in transgenic models of Alzheimer disease: implications for neural system disruption. J Neuropathol Exp Neurol 60:753–758
- Meyer-Luehmann M, Spires-Jones TL, Prada C et al (2008) Rapid appearance and local toxicity of amyloid-beta plaques in a mouse model of Alzheimer's disease. Nature 451:720–724
- Serrano-Pozo A, William CM, Ferrer I et al (2010) Beneficial effect of human anti-amyloid-beta active immunization on neurite morphology and tau pathology. Brain 133:1312–1327
- Spires TL, Meyer-Luehmann M, Stern EA et al (2005) Dendritic spine abnormalities in amyloid precursor protein transgenic mice demonstrated by gene transfer and intravital multiphoton microscopy. J Neurosci 25:7278–7287
- Susaki EA, Tainaka K, Perrin D et al (2014) Whole-brain imaging with single-cell resolution using chemical cocktails and computational analysis. Cell 157:726–739